



MR. Dental Supplies Limited

4 MANOR WAY, OLD WOKING, SURREY GU22 9JX. ENGLAND
Tel: +44(0)1483 773282 Fax: +44(0)1483 740548 Email: info@mr dental.co.uk Website: www. mr dental.co.uk

INSTRUCTIONS FOR USE

MEADWAY MODELLING COLDCURE POWDER:

Meadway Modelling Coldcure Powder is a Self-Curing Polymer for Various Applications in the Plastic, Dental, Film Industry.

Main uses are in the special effect departments of the film industry in the fabrication of Teeth, Fangs, Limbs, Joints etc.

MIXING:

The powder and liquid are mixed to a creamy consistency (1.3 grams powder to 1.0 ml of liquid)
Care should be taken when mixing to ensure that no air is trapped.

CURING:

It is not necessary to use a "pressure pot" to cure this material, and perfect results will be produced provided no air has been trapped during the mixing or building up stages, although a "pot" could be used if desired, to speed up the process.

The setting time in atmospheric conditions is 15-17 mins at 21 degrees C.

All other instructions are as per Meadway Rapid Repair and Meadway Premier Rapid Repair.

Liquid Used: Meadway Rapid Repair Coldcure Liquid, also known as Meadway Coldcure Modelling Liquid

WARNINGS AND CAUTIONS:

- *IN EXTREME RARE CASES MAY CAUSE SENSITISATION BY SKIN CONTACT*
- *USE ONLY BY TRAINED PROFESSIONALS*
- *IT IS THE RESPONSIBILITY OF THE TRAINED PROFESSIONAL TO ENSURE THAT THE PRODUCT IS SUITABLE FOR END USER*
- *STORE BETWEEN 16C TO 23C AWAY FROM DIRECT SUNLIGHT*
- *DISPOSABLE OF UNCURED LIQUID MUST BE CARRIED OUT IN ACCORDANCE WITH MATERIAL SAFETY DATA SHEET*
- *AND MIXED POWDER AND LIQUID CAUSES EXOTHERMIC REACTION*
- *IF THE DEVICE HAS BEEN PLACED IN THE ORAL CAVITY THEN IT MUST BE DISPOSED OF AS POTENTIALLY BIOLOGICAL HAZARD WASTE*
- *THE USER AND/OR THE PATIENT MUST REPORT ANY SERIOUS INCIDENT THAT HAS OCCURRED IN RELATION TO THE DEVICE TO THE MANUFACTURER AND THE COMPETENT AUTHORITY OF THE MEMBER STATE IN WHICH THE USER AND/OR PATIENT IS ESTABLISHED*



Issue No: ORIGINAL

Date: 21ST May 2021