



MR. Dental Supplies Limited

4 MANOR WAY, OLD WOKING, SURREY GU22 9JX. ENGLAND
Tel: +44(0)1483 773282 Fax: +44(0)1483 740548 Email: info@mr dental.co.uk Website: www. mr dental.co.uk

INSTRUCTIONS FOR USE

MEADWAY SPECIAL TRAY QUICK DOUGH POWDER- GREEN:

This is a cold cure material that must be used with its corresponding MEADWAY SPECIAL TRAY QUICK DOUGH LIQUID

1. Let the powder and liquid acquire the ambient temperature which should not be more than 23C
2. Pour 8.5 mls of the fluid into a mixing jar and add 16g of the powder to it.
3. Spatulate this mix for 30-45 seconds. The dough should be ready for application at this stage.
4. Transfer the dough immediately onto a polyethylene foil and spread it evenly to form a "plate".
5. Place this pre-formed "plate" on the plaster model and mould to prepare the functional tray. The dough should remain workable for 4-5 minutes, depending upon the ambient temperature.
6. Trim off the excess material along the periphery and use this to mould a handle which must be fixed on the tray while the dough is still workable.
7. Allow the tray to cure at the ambient temperature. The tray should be fully cured in a total of 8-12 minutes from the time of mixing, depending upon the ambient temperature.

WARNINGS AND CAUTIONS:

- IN EXTREME RARE CASES MAY CAUSE SENSITISATION BY SKIN CONTACT
- USE ONLY BY TRAINED PROFESSIONALS
- IT IS THE RESPONSIBILITY OF THE TRAINED PROFESSIONAL TO ENSURE THAT THE PRODUCT IS SUITABLE FOR END USER
- STORE BETWEEN 16C TO 23C AWAY FROM DIRECT SUNLIGHT
- DISPOSABLE OF UNCURED LIQUID MUST BE CARRIED OUT IN ACCORDANCE WITH MATERIAL SAFETY DATA SHEET
- AND MIXED POWDER AND LIQUID CAUSES EXOTHERMIC REACTION
- IF THE DEVICE HAS BEEN PLACED IN THE ORAL CAVITY THEN IT MUST BE DISPOSED OF AS POTENTIALLY BIOLOGICAL HAZARD WASTE
- THE USER AND/OR THE PATIENT MUST REPORT ANY SERIOUS INCIDENT THAT HAS OCCURRED IN RELATION TO THE DEVICE TO THE MANUFACTURER AND THE COMPETENT AUTHORITY OF THE MEMBER STATE IN WHICH THE USER AND/OR PATIENT IS ESTABLISHED

Issue No: ORIGINAL

Date: 21ST May 2021

